

Use the same cleaning technique as for a conventional manual toothbrush. Dentists recommend brushing teeth at least 3x a day for two minutes.



The practical change system enables you to replace the brush head at the latest after 3 months of use (recommended by dentists).



The TRISA Sonicpower is available with a rechargeable battery or with replaceable batteries.



The replaceable brush heads containing recycled plastic are an active contribution to the environment.







426423


MAW-vd2225-01 04/2022

TRISA AG | CH-6234 Triengen | Switzerland Phone +41 (0)41 935 35 35 | www.trisa.ch | bestellung@trisa.ch

For my brightest smile

TRISA Sonicpower







THE SONIC TOOTHBRUSH FOR EXTRA CLEAN TEETH

The TRISA Sonic power combines the comfort of use of a manual toothbrush with the advantages of a sonic toothbrush. Sonic vibrating bristles allow extremely effective teeth cleaning while protecting the gums in the best possible way. Developed and manufactured in Switzerland.



COMPLETE PROTECTION



The innovative «Complete Protection» brushing surface in combination with high frequency sonic movements cleans more than 30 per cent more efficiently.*



The scientifically developed «Complete Protection» brushing surface removes plaque with high efficiency and protects teeth as well as gums.





^{*} in comparison to conventional manual toothbrushes (ADA toothbrush)